

## PCR Thermal Cyclers Machine



**Product Categories:** [Lab Measuring Instruments](#)

**Product Page:** <https://www.labappara.com/product/2898/>

### Product Description

#### **PCR Thermal Cycler Machine - Scientific Equipment**

Adjustable pressure heat lid to prevent volatilization and dewing heat lid with pressure alarm device, to prevent damaging test tube by too much pressure. Temp, range 4°C to 99°C. Max. heating rate >2.5°C/s and Max. Cooling rate >2.5°C/s. Temp, display accuracy 0.1 °C. Max. cycle step 9. LCD Display graphical along with complete accessories and soft touch membrane keys operated.