

Heavy Duty Vacuum Bottles



Product Categories: [Bottles And Assesories](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page:

<https://www.labappara.com/product/42-heavy-duty-vaccum-bottles/>

Product Description

Heavy Duty Vacuum Bottles - Scientific Laboratory Equipment

The Reduction of convection provides the thermal insulation of thermos bottles.

Deep vacuum lowers the boiling point of liquids and promotes low temperature.

This autoclavable bottle, made of polypropylene, has thick walls, is rugged & can be used for repeated applications of vacuum.

Features:

Light weight

Crack resistance

Optimum strength

Smooth finish

Particulars	Packing
-------------	---------

2000 ml	1 Pc
---------	------