

Experiments to measure the Polar Pattern & gain characteristics of the antenna



Product Categories: [Electronics](#), [Engineering Equipment](#), [Microwave Lab Training Modules](#), [Microwave Lab Training Modules](#)

Product Page: <https://www.labappara.com/product/8033/>

Product Description

Experiments to measure the Polar Pattern & gain characteristics of the antenna

EXPERIMENTS :

- Pyramidal Horn
- Pick Up Horn
- Slotted W/G Horn
- Dielectric Horn
- Sectorial Horns
- Parabolic Dish

List Of Component:

- Klystron Power Supply
- Klystron Mount
- Klystron Tube
- Isolator
- Frequency Meter

Variable Attenuator

Radiation Pattern turn table

Standard gain horn

Detector mount

Cooling fan

VSWR meter solid state

Sectoral horn-E-plane

Sectoral horn-H-plane

Parabolic disc

Coaxial to W/G adopter

Twist

Moveable Short Precision

Stand

Bends

Tripod Stand

BNC Cable