

## Accelerated Curing Tank



**Product Categories:** [Pharmaceutical Equipments](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/accelerated-curing-tank/>

### Product Description

Accelerated Curing Tank - Scientific Equipment

Accelerated Curing Tank is used for concrete testing purposes.

Its inner chamber is made of stainless steel and outer body is made of mild steel which is duly powder coated.

Suitable to accommodate 6 cube moulds of 150 mm.

The unit is provided with digital temperature controller cum indicator and water circulation pump.

Temperature range 5° above ambient to 100°C.

Calibration certificate is provided for digital controller.

Special Features :

Temperature Controller Cum Indicator

Water Circulation Pump

Calibration Certificate is Provided For Digital Controller