


# Adapter Borosilicate Glass



**Product Categories:** [Adapters Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/adapter-borosilicate-glass/>

## Product Description

Adapter Borosilicate Glass - Laboratory Glassware Equipment

Adapter borosilicate glass with interchangeable joints are manufactured from finest quality borosilicate glass which gives maximum resistance to thermal and mechanical stress and chemical attack. Multiple adapters with three necks, two neck parallel and one at 45°. Adapter is made of borosilicate glass and joints are interchangeable with other apparatus.

Socket Size	Cone
14/23	14/23
19/26	19/26
19/26	24/29
24/29	24/29
19/26	34/35
29/32	29/32