

Amber Colour Burettes



Product Categories: [Burettes Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [borosilicate glass](#), [laboratory glassware](#)

Product Page: <https://www.labappara.com/product/amber-colour-burettes/>

Product Description

Amber Colour Burettes- Laboratory Glassware Equipments

Automatic zero, with Glass Key Stop Cock mounted on reservoir, accuracy as per 'Class A' with rubber bellow.

Available In Capacity(ml)- 10 ml, 25 ml, 50 ml, 100 ml

Available In Sub-div (ml)- 0.05, 0.1, 0.1, 0.2