


Amber Color Pipettes



Product Categories: [Laboratory Glassware](#), [Pipettes](#) , [Volumetric Flasks And Measuring Cylinder](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/amber-colour-pipettes/>

Product Description

Amber Colour Pipettes - Laboratory Glassware Equipment

A pipette is a laboratory tool commonly used in laboratory glassware a measured volume of liquid, often as a media dispenser

Pipettes come in several designs for various purposes with differing levels of accuracy, from single piece glass pipettes to more complex adjustable or laboratory pipettes

Features:

Durable finish

Best quality

Premium grade materials are used

Specifications:

Vol.: 1ml

Material: Glass

"Labappara", Accuracy as per I.S. 4162:

Cat. No.	Capacity	Sub Div.	Tolerance
ml.	ml.	ml.	
39A/1	0.1	0.01	0.006
39A/2	0.2	0.01	0.006
39A/3	0.5	0.02	0.006
39A/4	1.0	0.01	0.006
39A/5	1.0	0.1	0.006
39A/6	2.0	0.02	0.01
39A/7	2.0	0.1	0.01
39A/8	5.0	0.1	0.03
39A/9	5.0	0.05	0.05
39A/10	10.0	0.1	0.05
39A/11	25.0	0.2	0.1