

Ampoule Filling & Sealing Device



Product Categories: [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/ampoule-filling-sealing-device/>

Product Description

Ampoule Filling & Sealing Device - Laboratory Equipment

Electrically operated with piston of syringe having adjustment for quantity of liquid to be filled.

Completed with burner for sealing.

Features:

Unique designs

Exclusiveness

Appealing look

Complete motorized capacity adjustments

Longer service life

Power efficiency