

## Analog Lab Trainer (With Breadboard)



**Product Categories:** [Bread Board Trainers](#), [Bread Board Trainers](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/analog-lab-trainer-with-breadboard/>

### Product Description

Analog Lab Trainer (With Breadboard)

#### **Specification:**

Trainer comprises of Power supplies inside the trainer

Power supplies on boards 0-50V, 250 mA and 0-30V,250mA Variable  
15V/150mA and 5V/150mA Fixed Supplies provided on board.

4 Nos of Variable resistor provided on the board.

One 20 Pin Vero board on mounted on the board.

One Bread Board mounted on board.

Ac power supply 15-0-15V for rectifier applications

Various components on the trainer set

3 Pairs of connected 2mm sockets

4 Rectifier diodes

1 OA79 germanium diode

3 different values of Zener

N different values resistances

N different values capacitors

Inductors 1mh, 10mh

Thyristor Bank

SL100 and SK100 pair

1 LDR

1 Thermistor

2 LED's

1 Relay 12V

2 Nos of Digital Voltmeter and Ammeter

One Speaker 8E 1/4W

1 Audio Transformer