

Astable Multivibrator (IC 555)



Product Categories: [Electronics](#), [Engineering Equipment](#), [Multivibrator & Oscillator](#), [Multivibrator & Oscillator](#)

Product Page: <https://www.labappara.com/product/astable-multivibrator-ic-555/>

Product Description

Astable Multivibrator (IC 555)

555 Timer as Astable Multivibrator Apparatus has been designed to study the Wave Shape and Frequency produced by Astable Multivibrator.

The instrument comprises of the following built in parts:-

IC 555 placed inside the cabinet & connections brought out at sockets.

DC regulated power supply of 5V.

Various resistances & capacitors connected inside the cabinet & connections brought out at the sockets.