

## Audio Oscillator Solid State



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:** <https://www.labappara.com/product/audio-oscillator-solid-state/>

### Product Description

Audio Oscillator Solid State - Electronics Laboratory Instrument

#### **Features:**

Accurate output

Easy operation

Superb strength

#### **Other Details:**

(Sine & Square Wave)

Audio oscillator solid state Frequency Range - 1 Hz to 100 KHz in 5 decade steps

Output Voltage 0 - 10 V p-p Continuously Variable

Selector Switch for Amplitude 0-1 / 0-10 V is provided

Selector Switch for Sine & Square Wave is provided

Output Impedance - 50 ohms