

## Auto Conductivity Meter



**Product Categories:** [Analytical Lab Instruments](#), [Conductivity Meters](#), [Laboratory Equipments](#)

**Product Tags:** [Analytical Lab Equipment](#)

**Product Page:** <https://www.labappara.com/product/auto-conductivity-meter-2/>

### Product Description

Auto Conductivity Meter -Analytical Lab Equipment

Alphanumeric LCD Display with facility for automatic temperature compensation and digital.

Cell constant adjustment along with conductivity cell & temperature probe.

Soft touch membrane keys.

- a) Conductivity Cell Stand
- b) Conductivity Cell (Cell constant 0.1)
- c) Conductivity Cell (Cell constant 0.5)
- d) Conductivity Cell (Cell constant 1.0)

e) Standard solution for conductivity 1408/vS or 1000//S .

f) Temperature Prob

Auto Digital Conductivity Meter (Auto Ranging) 3½ digit LED display with facility for temperature compensation and digital cell constant adjustment along with conductivity cell & other accessories (Table Model)

Range	5 Ranges 0 to 1000m Mhos/cm
Display	3½ Digit LED Display
Accuracy	± 2% FS ± 1 Digit
Resolution	0.1 micro mhos/cm
Comp.	0° to 50° C manual
Cell Constant	Adjustable on display
Cell	Platinum Dip
Accessories	Conductivity Cell, Dust Cover, Manual.