

Auto Conductivity-Tds-Temp Meter



Product Categories: [Analytical Lab Instruments](#), [Conductivity Meters](#), [Laboratory Equipments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/auto-conductivity-tds-tempmeter/>

Product Description

Auto Conductivity-Tds-Temperature Meter - Analytical Laboratory Equipment

Alphanumeric LCD Display with facility for automatic temperature compensation and digital cell constant adjustment along with conductivity cell, IDS Cell & temperature probe.

Soft touch membrane keys.

Function

Cond. & Temperature , TDS & Temperature, SAL & Temperature

Simultaneously

Alphanumeric LCD Display with facility for automatic Temperature compensation and digital cell constant adjustment along with conductivity cell, IDS Cell & temperature probe.

Soft touch membrane keys.

Contains1 : Conductivity/TDS cell, temperature probe, instruction manual

Display Text : 16 x 2 alphanumeric LCD display

Frequency : 50 Hz

Range : 0.20, 0-200 us, 0-2, 0-20, 0-200 ms

Resolution : 0.01 us

Storage : Data Sample: up to 100 Samples

Temperature Range : auto/manual: 0 to 50??C

Voltage : 230 V +10% AC