

Auto Potentiometer



Product Categories: [Analytical Lab Instruments](#)

Product Page: <https://www.labappara.com/product/auto-potentiometer/>

Product Description

Auto Potentiometer - Analytical Laboratory Equipment

Redox Potential Measurement.

Soft touch membrane

8 keys display 4 digit LED range 0 to ± 1999 mV.

Power supply 230V \pm 10% supplied with Glass, References (Calomel), Silver & Platinum Electrodes.

Range

pH

0 to 200

mV

0 to 2000

Accuracy

pH

mV

$\pm 1 \text{ mV} \pm 1 \text{ digit}$

Operating temperature

10° c to 45° c

Slope correction

80% to 120%

Ambient Temperature

0° C to 50° C

Display

3½ digit LED display

Power

230 V \pm 10% 50Hz AC

Dimensions

275x175x85 mm

Accessories

1) Magnetic stirrer2), Silver and Platinum Electrodes

3) Dust cover

4) Instruction manual.

Platinum electrode : 1 no.

Reference electrode : 1 no.

Glass electrode: 1 no.

Silver electrode : 1 no.

Dust cover: 1 no.

Electrode stand: 1 no.