

## Autoclave (Electrical ,Cooker Type ,S.S)



**Product Categories:** [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/autoclave-electrical-cooker-type-s-s/>

### Product Description

Autoclave (Electrical ,Cooker Type ,S.S) - Laboratory Equipment

Made of Stainless Steel. Also available in Electric / Non Electric (Fuel Based). Also available with Semi-Automatic (Timer Based). All sizes are approximate.

Used for Sterilization of Surgical Instruments and Dressings.

Seamless Construction (Pointless) from Extra Thick Aluminum Sheet/S.S. Sheet for better safety

Good quality synthetic rubber Gasket for Longer and better performance and Hygiene.

Weight Type Safety Valve and Pressure Regulating Device for better safety and functioning.

Co loured ceded Pressure gauge (0-30 bis/in<sup>2</sup>) to indicate the safe, sterilization and danger zones.

Highly strong and durable vacuum and steam resale valves.

Handles and valve knobs made of very durable, heat resistant and fine quality Bakelite.

Heating element made of heavy copper pipe for longer life.

Well designed, powder coated outer stand to provide better stability and attractive looks.

Well Packed in Individual boxes.

**Size 12" x 12"**