

Autoclave (Gas)



Product Categories: [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/autoclave-gas/>

Product Description

Autoclave (Gas) - Laboratory Equipment

Portable autoclave is made up of three walls and the three walls are made of high grade stainless steel sheet.

Suitable for natural/bottled gas.

The container is made of seamless aluminum body and lid.

A rubber gasket makes the lid and body airtight with help of tightening bolts and nuts.

Lid is fitted with pressure gauge and double safety valve.

Thermally insulated and vulcanized poly chloro-prone (Synthetic Rubber) gasket.

The construction of the sterilizer is ensured an easy and safe handling.

Handles holding lugs 7 nobs afford safe handling in-spite of high temperature.