

AUTOCLAVE (NON ELECTRICAL)



Product Categories: [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/autoclave-non-electrical/>

Product Description

AUTOCLAVE (NON ELECTRICAL) - Laboratory Equipment

Suitable to sterilize bandages & needles. Gas or Stove heated. All thick aluminum body with pressure gauge and safety valve. These autoclaves are made of 10 gauge aluminum. The drum is deep drawn & the ring on top is riveted. The lid is made by casting. Locking of lid by 6 nuts with rubber gasket seal in between (not metal to metal seal). Two safety valves, pressure gauge, flexible metal rubbing to exhaust air. Includes outer stand, inner aluminum container, inner tripod stand.

Size 12" x 12"

Size 12" x 15"