

Autoclave (Portable Single Drum)



Product Categories: [Autoclave & Sterilizer](#), [Autoclaves](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page:

<https://www.labappara.com/product/autoclave-portable-single-drum/>

Product Description

Autoclave (Portable Single Drum) - Scientific Equipment

Portable single drum Autoclave, highly efficient apparatus designed to meet the requirements of Labs.

Very simple to operate and economical to use.

The container is seamless made of thick stainless steel.

The unit is fitted with radial locking system with paddle lifting device, made of MS chrome plated.

Provided with built- in safety valve, pressure gauge, pressure release valve and water level indicator.

The pressure is adjustable from 5 psi to 20 psi with automatic pressure control switch, which controls the pressure with an accuracy of +/- 1 psi.

Supplied complete with stainless steel basket, chord and plug.

