

## Auxanometer



**Product Categories:** [Plant Physiology Apparatus](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/auxonometer/>

### Product Description

#### Auxanometer - Laboratory Equipment

An auxanometer is an apparatus for measuring increase or rate of growth in plants. In case of an arc-auxanometer there is a wire fixed with the plant apex on one end and a dead-weight on the other. It passes over a pulley which has a pointer attached to it. When the plant's height increases, the pulley rotates and the pointer moves on a circular scale to directly give the magnitude of growth.

Auxanometer is with new clock movement having start and stop system.

An ideal instrument for recording growth of plants.

Auxanometer is complete with aluminum pulley, plate, brass clamp and recorder.

And within six rods fitted on heavy cast iron base.