

## B. H Curve Apparatus



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:** <https://www.labappara.com/product/b-h-curve-apparatus/>

### Product Description

B. H Curve Apparatus - Electronics Equipment

B.H curve to trace hysteresis curves on C.R.O.

Terminals are provide for x & y plates.

B.H using resistance , capacities & inductance step down transformer 2-12 volt AC.

#### **Features:**

Instrument comprises of the following:

Output selectable using band switch.

Circuit diagram for BH curve printed.

Components connected behind the front panel.

Input & output connections brought out at sockets.

One ON/OFF Switch with jewel light is provided on the front panel.

Made of Heavy duty metal box construction.

Front Panel made of high quality Bakelite sheet.

Strong construction

Easy operations

High strength

**Accessories Supplied:**

Required number of Patch Chords.

Instruction Manual.