

Basic Mobile Trainer



Product Categories: [Audio Visual Lab Demonstrators & Trainers](#), [Audio Visual Lab Demonstrators & Trainers](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/basic-mobile-trainer/>

Product Description

Basic Mobile Trainer

FEATURES :

It has hardware block for doing various experiments

Mic input with conditioning circuit

Keyboard matrix (4x4) for scanning the key codes

Audio output on Headset

SIM socket

Buzzer

With buzzer provides for incoming call

External antenna for capturing the incoming signal

In-built power supply

SPECIFICATIONS :

Cellular System : EGSM/GSM 900

Rx Frequency band : EGSM 925..960MHz GSM 900 935 .960MHz

Tx Frequency band : EGSM 880 ..8790MHz GSM 90 890 915 MHz

Output power : +5 +33dBm/3.2 mW ..2W

Channel Spacing : 200KHz

Antenna : Loop type, 50W

Display : 84x48 pixels

On board sections : Antenna, Keypad, SIM, Charging Circuit Clock, User Interface such as Buzzer, Vibrator, LEDs.

Features : Screen saver, Ring tones, SMS etc.

Accessories included : Battery, Mains cord, Manual, Hands Free Kit CD containing Mobile Phone working

Power requirement : 220V \pm 10% 50Hz

Power consumption : 3.6 Watts (Approx.)

Fuse : 1.5 amps

Housed in attractive cabinet