

Bath (Sand)



Product Categories: [Laboratory Equipments](#), [Water Bath& Stills](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/bath-sand/>

Product Description

Bath (Sand) - Scientific Laboratory Equipment

Shallow form, made of thick G.I. sheet. Double wall construction, Inner .made of stainless steel. Outside cement sheet fixed on M.S. angles. The bath is heated by open type heaters controlled by energy regulators.

Features:

Compact size

Robust construction

Excellent functioning