

Biological Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/biological-microscope/>

Product Description

Biological Microscope - Biology Laboratory Equipment

The Biological Microscope are special designed for education of students for study in Biology, Microbiology, Anatomy departments of Universities, Medical Colleges and Hospitals. Provided with Achromatic Objectives 4 x/5 x, 10 x, 40 x & 100 x Oil and eyepieces 10 x. The magnification range is 20 x to 1500 x with different objective and eyepieces. It is available in three models that are Monocular Biological Microscope, Binocular Biological Microscope and Binocular Biological Microscope with digital camera.

Features:

Foldable eye guards

Precisely designed

Reliable usage

Specification: Stand: Very light weight stand with single focus control and stage

clips and small reflector

Optics: Eye piece 10 X or 15 X, objective 4 X or 10 X or 20 X