

Blender



Product Categories: [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/blender/>

Product Description

Blender - Laboratory Equipment

Blender is electrically operated

Speed can be controlled with speed control device

With two stainless steel jars of different sizes.

Motor power (w) : 550 W.

Revolution : 18000 RPM.

Power requirement : 230 V, 50 Hz.

Unbreakable and dishwasher safe.

It has easily removable container cover with removable clear cap for adding ingredients and for pressure reduction during hot food blending.