

# Blood Donor Couch (Electrical Blood Couch)



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/blooddonor-couch-electrical-blood-couch/>

## Product Description

Blood Donor Couch (Electrical Blood Couch) - Hospital Equipment

Facility for blood collection from both sides. Variable positioning for either arm with comfortably wide arm- rests.

Two geared Actuator.

Based on homo dynamic principles. (Donor's head can be lowered immediately and legs can be lofted above vital organs)

4 Heavy Duty Lockable castors for easy mobility.

Comfortable chair type with soft padding for cushioning and rexcin cover.