

## Bottle Borosilicate Glass



**Product Categories:** [Bottles Borosilicate Glass](#), [Laboratory Equipments](#), [Quartz Ware](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/bottle-borosilicate-glass/>

### Product Description

Bottle Borosilicate Glass - Laboratory Glassware Equipment

Borosilicate glass is a type of glass and boron trioxide as the main glass-forming constituents.

Borosilicate glasses are known for having very low coefficients of thermal expansion.

Borosilicate glasses them resistant to thermal shock, more so than any other common glass.

#### **Features:**

Hygienic and clean

Easy to dispose