

## Bottle / Liquid Filling Machine



**Product Categories:** [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/bottle-liquid-filling-machine/>

### Product Description

Bottle / Liquid Filling Machine - Laboratory Equipment

A sturdy and heavy unit made of sheet mounted on a heavy angle iron duly powder coated.

Capacity of output of water is adjustable through cam system from 50 to 300 ml.

#### **Technical Specifications:**

Filling capacity variations: 50 to 300 ml

Syringe options: Single/double

Stainless steel body

Semi automatic

Electric operated

#### **Features:**

Rust proof body

Easy to operate

Long service life