

Brain Model



Product Categories: [Anatomical Models](#), [Biology Lab](#), [Educational equipments](#), [Educational Models](#), [Organ Models](#)

Product Tags: [Scientific Educational Model](#)

Product Page: <https://www.labappara.com/product/brain-model/>

Product Description

Brain Model - Scientific Educational Model Brain Model-Brain features of human brain as a whole, as well as the relations of between the component portions. Dissectible into 3 parts. Natural size. Made of unbreakable high quality PVC and hand painted with non-toxic colors in a scientifically correct manner. Separates into 2 parts, median section. Mounted on a plastic base, length approx. 200 mm, width approx. 165 mm, thickness approx. 150 mm. It is supplied with detailed instruction cards and conforms to the main EU safety standards.

Features:

Crack resistance

Impeccable finish

Light weight

Size: 18.4x14x13.5 cm.