

## Buchner Funnel



**Product Categories:** [Funnels Borosilicate Glass](#), [Laboratory Glassware](#), [Sintered Ware Borosilicate Glass](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/buchner-funnel-2/>

### Product Description

Buchner Funnel - Laboratory Glassware Equipment

Buchner Funnel with Sintered Disc of Porosity G-1 or G-2 or G-3 or G-4, plain stem. Moulded in Polypropylene, these two piece funnels are light weight, easy to clean & autoclavable. The top portion can be snap fitted with the lower one forming a vacuum seal which does not require greasing. The base of the top portion has a perforated filter plate.

#### Features:

Rugged in construction

Light in weight

Reliable performance

Available in Capacity(ml) - 35 ml, 80 ml, 200 ml, 500 ml, 1000 ml.

Available in disc dia (mm)- 30 mm, 40 mm, 65 mm, 90 mm, 120 mm.