

# Burette Brush



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/burette-brush/>

## Product Description

### Burette Brush

To hold one burette and one funnel.

Polished hard wood with cork lined jaws & plastic screws.

Base fitted with four rubber feet for stability.

Overall height 490mm,  
base size 330x110mm.

Reliable operation

Low maintenance

Easy to operate