

## Burette Clamp (Fisher Type)



**Product Categories:** [Chemistry Equipments](#), [Lab Supplies](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/burette-clamp-fisher-type/>

### Product Description

Burette Clamp (Fisher Type) - Chemistry Equipment

**Burette clamp** is extensively used for igniting bunsen burners to gas supply.

#### **Feature -**

Burette Clamp Fisher Type with improved design and superb finish, made of thick mild steel sheet, die pressed, symmetrical non breakable with two jaws having 'V' grooves in which Burette rest and two spring controlled straight jaws with even & precise pressure hold the Burette safely.

Comes with Die pressed mild steel Boss Head with 'T' type screw which accommodate up to 18 mm dia rod.

Comes with special Powder Coated superb finish.

Power coated die pressed non breakage.

For single burette.

For double burette.

Complete with boss head die pressed iron packed.

Brass Oxidized with boss head.