

Burettes with Screw Type PTFE Stopcock



Product Categories: [Burettes Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [biology lab equipment](#), [chemical glassware](#), [chemistry lab instruments](#), [Industrial Glassware](#), [lab glassware](#), [scientific glassware](#)

Product Page:

<https://www.labappara.com/product/burettes-screw-type-ptfe-stopcock/>

Product Description

Burettes with Screw Type PTFE Stopcock - Laboratory Glassware Equipment

Made of Borosilicate glass.

Graduations and inscriptions in blue enamel

Schellbach burette tube with vertical white stripe and central blue ribbon

Exact position of the meniscus is clearly defined where the blue ribbon splits into a wide and narrow band

With funnel top to make filling easier and safer

Each burette has a +30 s waiting time before final reading can be taken

Grease-free interchangeable PTFE key

Features:

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

Auto-zero on Reservoir Plain with Stop-cock

| Capacity | Sub. Div. |
|----------|-----------|
|----------|-----------|

ml

| | |
|----|-----|
| 10 | 0.1 |
|----|-----|

| | |
|----|-----|
| 25 | 0.1 |
|----|-----|

| | |
|----|-----|
| 25 | 0.1 |
|----|-----|

| | |
|----|-----|
| 25 | 0.1 |
|----|-----|

| | |
|----|-----|
| 50 | 0.1 |
|----|-----|

| | |
|----|-----|
| 50 | 0.1 |
|----|-----|

| | |
|-----|-----|
| 100 | 0.2 |
|-----|-----|