

Cage Metabolic



Product Categories: [Laboratory Equipments](#), [Pharmaceutical Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/cage-metabolic/>

Product Description

Cage Metabolic - Scientific Equipment

Metabolic Cage : This metabolic cage is used for the efficient collection of urine and feces from a single mouse. The smaller chamber volume and collection tube design minimize evaporation. The funnel and cone design, with a steeper angle, separates urine and feces into separate collection tubes. Non-wetting funnel and cone are made of polymethylpentene (PMP) plastic to ensure immediate and complete separation of urine and feces. Cage is self-standing on a bench-top or shelf rack.

Specification:

Floor Area	:	200 cm ² (31 in. ²)
Chamber Volume	:	2600 cm ³ (159 in. ³)

Dimensions Cage : 32 L x 25 W x 36.5 H cm

Vertical Support Rack : 126 L x 50 W x 187 H cm

Horizontal Support Rack : 178 L x 63 W x 129 H cm

Features:

Reliable results of experiments

Practicality of use

Flexibility

Durability

Space saving