

Cam Analysis Apparatus



Product Categories: [Engineering Equipment](#), [Theory Of Machine Lab](#)

Product Page: <https://www.labappara.com/product/cam-analysis-apparatus-2/>

Product Description

Cam Analysis Apparatus

Description:

With the help of our dexterous professionals we have come up with Cam Analysis Apparatus. These are provided with a DC Motor and an extended shaft that carries a cam. The procuring agents of our firm ensure that the raw material for the product is procured from certified vendors only. To check the compliance of the product with international quality standards, these are stringently tested by our quality analyzers, post production.

Features:

Easy operations

Dimensional accuracy

Energy efficient

Accurate measurement

Durability

Precisely designed

Specifications:

Cam -Eccentric, tangent and circular ARC type

Follower- mushroom, flat faced and roller type

Cams and followers are hardened to reduce wear of the surfaces

Variable speed motor coupled to camshaft of suitable range

A dial gauge to note the follower displacement

A technical manual accompanies the equipment