

Cascaded Two Port Network



Product Categories: [Electronics](#), [Engineering Equipment](#), [Filter Circuits Lab](#), [Filter Circuits Lab](#)

Product Page: <https://www.labappara.com/product/cascaded-two-port-network/>

Product Description

Cascaded Two Port Network

Specification :

One continuously variable, overload and short circuit protected dc regulated power supply of 0- 10v. Two dual range meters are mounted on the front panel to measure the voltage & current. Circuits for two networks are given on the front panel (resistances in each arm are of 1kw.)