

Casseroles



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/casseroles/>

Product Description

Casseroles - Scientific Laboratory Equipment

Porcelain Casseroles, Basins with handles. Rugged, with pouring lip and handle. Glazed inside out except on rim and bottom. Flat stable bottom. All dimensions are approximate.

Diameter :- 50 mm, 75 mm, 85 mm, 100 mm, 150 mm, 175 mm

To prolong the life of porcelain products:

Follow gradual heating/cooling rates

Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required

Avoid contact of heated ware with a cold surface

Carefully inspect lab ware prior to each use and do not use any product that appears defective