

## Centrifuge (Economical)



**Product Categories:** [Centrifuges Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-economical/>

### Product Description

Centrifuge (Economical) - Scientific Laboratory Equipment

Economical unit for small capacity centrifuging requirements. The standard unit is fitted with resiliently mounted motor with 3500 RPM for vibration free performance, a built in 5 step speed regulator and one non removable angle motor. Complete with tube carriers.

**Features:**

- Compact size
- Sturdy design
- Flawless finish
- Long working life
- High functionality

**Available with following angle head rotors -**

- 15 ml x 4 tubes
- 15 ml x 6 tubes
- 15 ml x 8 tubes