

Centrifuge (Hand 2 tubes)



Product Categories: [Centrifuges Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/centrifuge-hand-2-tubes/>

Product Description

Centrifuge (Hand 2 tubes) - Scientific Laboratory Equipment

A single speed hand driven machine consisting of a die casted parts and a steel pinion and worm drive. Screw clamps the machine to the bench. Well Finished , with removable crank handle. Suitable for 15 ml glass tubes.

This hand driven centrifuge comes with a table-mounting clamp. Features a painted metal body, and metal head.

Comes complete with 4 each aluminum shields.

Easily dis-assembles for storage in specially designed foam box.