

## Centrifuge (Laboratory)



**Product Categories:** [Centrifuges Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-laboratory/>

### Product Description

Centrifuge (Laboratory) - Scientific Laboratory Equipment

Microprocessor based. Square M S body duly powder coated. Double walled light weight ABS lid. Fitted with microprocessor based 2 lines 16 characters LCD panel for 0-59 minutes countdown timer, digital rpm meter and programmable speed controller.

For Clinical & Pathological Laboratories. Suitable for infertility labs where low speed is recommended.

Brush less AC Motor

Sea thru acrylic lid.

Digital Timer 0-99 minutes.

Maximum speed 3000 r.p.m. without load and 2850 rpm on load

With 8 x 15ml Swing out head

With 8 reduction adopter to accept test tubes of size 12 x 100 mm.

Laboratory Equipments for use in Schools, Colleges, Laboratories, Industries and

Research Institutes.

**Capacity:**

4x15 ml tubes.

6x15 ml tubes.

8x15 ml tubes.

12x15 ml tubes.