

Centrifuge (Mini)



Product Categories: [Centrifuges Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/centrifuge-mini/>

Product Description

Centrifuge (Mini) - Scientific Laboratory Equipment

Maximum Speed 10000 rpm Spin-Top with 2 heads 6 x 1.5 ml and 12 x 0.5 ml tube. Microprocessor based digital timer, digital rpm and speed control from 1000 rpm to 10000 rpm in 10 steps. Fitted with imported motor from USA without carbon.

Mini Centrifuge Digital

Micro processor controlled, Micro-Centrifuge machine (Mini Centrifuge)

Speed range from 1000 to 10000 r.p.m. in step of 1000.

Technical Specifications:-

Fitted with maintenance free motor (Brushless without carbon) Complete body made of ABS plastic injection moulded.

Fitted with 0-15 minute digital timer.

Supply completed with 2 standard rotors as under :

- (i) Angle Rotor 6 micro centrifuge of 1.5ml capacity
- (ii) Angle rotor 12 micro centrifuge tubes of 0.5 ml capacity