

## Centrifuge Tube Box



**Product Categories:** [Centrifuges Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-tube-box/>

### Product Description

#### Centrifuge Tube Box

Radical Centrifuge tube box are molded in polypropylene. These boxes provide excellent substitute for conventional card-box or foam rack. The Transparent top gives a clear view of the tubes placed in the box which are easily identifiable because of the indexed grid. These polypropylene centrifuge boxes are a more durable alternative to foam racks. Boxes are auto cleanable.

#### **Features:**

Crystal clear lid

Maintain integrity under harsh freezer condition

These racks can withstand a temperature of -90° C to 121° C.

The distinctive green color on the 15 ml and blue color on the 50 ml makes these racks easy to distinguish.

All grids are properly indexed with clear lid that provides easy tube identification.

Can be autoclaved.

Capacity:- 15 ml

Capacity:- 50 ml