

Centrifuge Tube Round Bottom



Product Categories: [Centrifuge Bottles And Tubes](#), [Laboratory Plastic Ware](#)

Product Tags: [Labaoratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/centrifuge-tube-round-bottom/>

Product Description

Centrifuge Tube Round Bottom - Scientific Laboratory Equipment

Conical Bottom Centrifuge Tube: Simplify measurements using Conical-Bottom Centrifuge Tubes, which feature molded-in graduations. Translucent polypropylene co-polymer (PPCO) has excellent chemical resistance. These tubes can be autoclaved, however, repeated autoclaving in PC tubes reduces their mechanical strength.

Particulars	Packing
16 ml. Cent Tube PP	50 Pcs
50 ml. Cent Tube PP	50 Pcs
16 ml. Cent Tube PC	25 Pcs
50 ml. Cent Tube PC	25 Pcs
Cap for 16 ml. Cent Tube	100 Pcs
Cap for 50 ml. Cent Tube	100 Pcs