

## Centrifuge Tube



**Product Categories:** [Test Tubes / Culture Tubes](#), [Test Tubes / Culture Tubes](#)  
[Borocilicate Glass](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-tube/>

### Product Description

Centrifuge Tube - Laboratory Glassware Equipment

"Labappara", Conical Bottom, graduated

High quality, crystal-clear polypropylene centrifuge tubes. Centrifuge tubes can be safely centrifuged to approximately  $3,200 \times g$  (50 m L size) and  $3,500 \times g$  (15 m L size) depending on the specific gravity of the solution being centrifuged and the type of rotor.

#### **Features:**

Perfect finish

Thermal resistance

Resistant to crack

Available in capacity - 5 ml, 10 ml, 15 ml, 20 ml, 25 ml.