

## Centrifuge Tubes, Oil



**Product Categories:** [Laboratory Glassware](#), [Test Tubes / Culture Tubes](#), [Test Tubes / Culture Tubes Borocilicate Glass](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-tubes-oil/>

### Product Description

Centrifuge Tubes, Oil - Laboratory Glassware Equipment

High quality, crystal-clear polypropylene centrifuge tubes. Centrifuge tubes can be safely centrifuged to approximately 3,200 × g (50 m L size) and 3,500 × g (15 m L size) depending on the specific gravity of the solution being centrifuged and the type of rotor.

#### **Features:**

Perfect finish

Thermal resistance

Resistant to crack

Fully graduated 0 to 0.5 ml.in 0.05 ml.

Div., 0.5 to 2 ml. in 0.1 ml.

Div., 2 ml. to 3 ml. In 0.2 ml.

Div., 3 ml. to 5 ml. in 0.5 ml.

Div., 5 ml. to 10 ml. in 1 ml.

Div., 10 ml. to 25 ml. in 5 ml.

Div., graduations are marked at 50, 75 and 100 ml.