

Centrifuge Tubes (Round bottom)



Product Categories: [Centrifuges Equipment](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/centrifuge-tubes-round-bottom/>

Product Description

Centrifuge Tubes (Round bottom) - Scientific Laboratory Equipment

Radical Centrifuge Tubes with round bottom are molded in Polypropylene & Poly carbonate These tubes have excellent clarity and are strong enough to be used up to 50,000 x G. These tubes can be autoclaved; however, repeated autoclaving in PC tubes reduces their mechanical strength.

Avail from us, a flawless variety of Round Bottom Centrifuge Tube, which is popular for its better quality and precise engineering.

Features:

Centrifuge tube with white enamel graduations

Graduated 0 – 3 ml in 0.1 ml divisions and 3 – 10 ml in 0.5 ml divisions

Can be centrifuged up to 3000 RCF values

Centrifuge Tube Round Bottom 15 ml (P.P)

Centrifuge Tube Round Bottom 50 ml (P.P)