

# Centrifuge (Universal)



**Product Categories:** [Centrifuges Equipment](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/centrifuge-universal/>

## Product Description

Centrifuge (Universal) - Scientific Laboratory Equipment

This Centrifuge Compact model for routine work in laboratories. The special features include easy lid lock, sleepless speed regulator and timer.

**Specification:**

Lab Centrifuge is Electrically operated and Works on 230/50Hz volts

Speed regulation is possible using rotational knob

Attractive enamel bacteria resistant powder coated finishing of Lab Centrifuge

Meditech Lab Centrifuge works on brushless motor with 5000 R.P.M

Rotor : swing out head and angle design

Tubes : 4 tubes, 8 tubes and 12 tubes configuration.

**Features:**

Excellent performance

High durability

Longer functional life

It is extensively used in bio technology, hospitals, blood banks, educational and research institutions, pharmaceutical and medical laboratories for centrifugal of temperature sensitive material.

**Available with following rotors :-**

15 ml x 4 tubes, 5000 RPM, Swing out head

15 ml x 6 tubes, 5000 RPM, Swing out head

15 ml x 8 tubes, 5000 RPM, Swing out head

15 ml x 12 tubes, 4200 RPM, Angle head

50 ml x 4 tubes, 2800 RPM, Angle head

50 ml x 6 tubes, 2400 RPM, Angle head

100 ml x 4 tubes, 2500 RPM, Angle head