

China Dish (Round form)



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/china-dish-round-form/>

Product Description

China Dish (Round form) - Scientific Laboratory Equipment

This porcelain evaporating dish is thermally-stable and glazed inside for easy cleaning. Heat chemical solutions in a porcelain evaporating dish, driving off water and leaving chemical solute behind

Features:

Crack resistance

Easy to clean

Can withstand high temperature

| Diameter | capacity |
|----------|----------|
| 25 mm | 15 ml |
| 35 mm | 20 ml |
| 50 mm | 40 ml |
| 60 mm | 50 ml |
| 75 mm | 100 ml |
| 85 mm | 125 ml |

| | |
|--------|----------|
| 100 mm | 175 ml |
| 110 mm | 275 ml |
| 125 mm | 450 ml |
| 150 mm | 650 ml |
| 200 mm | 1000 ml |
| 250 mm | 2000 ml |
| 300 mm | 3000 ml |
| 350 mm | 6500 ml |
| 400 mm | 12000 ml |
| 450 mm | 16000 ml |
| 500 mm | 18000 ml |