

Chromatography Assembly Semi Micro



Product Categories: [Chromatography Vials Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipments](#)

Product Page:

<https://www.labappara.com/product/chromatography-assembly-semi-micro/>

Product Description

Chromatography Assembly Semi Micro - Laboratory Glassware Equipment

Chromatography Assembly Semi Micro consisting of funnel Reservoir 50 ml. capacity with 14/23 cone and socket,

Column tube with tap, 12 mm. dia 10 cm.

Effective length with 14/23 socket, 19/26 cone and fused sintered disc.

Flask 100 ml. with 19/26 socket tamping rod 40 cm. Long.

Features:

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

Auto-zero on Reservoir Plain with Stop-cock