

## Chromatography Assembly Semi Micro



**Product Categories:** [Chromatography Vials Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipments](#)

**Product Page:**

<https://www.labappara.com/product/chromatography-assembly-semi-micro/>

### Product Description

Chromatography Assembly Semi Micro - Laboratory Glassware Equipment

Chromatography Assembly Semi Micro consisting of funnel Reservoir 50 ml. capacity with 14/23 cone and socket,

Column tube with tap, 12 mm. dia 10 cm.

Effective length with 14/23 socket, 19/26 cone and fused sintered disc.

Flask 100 ml. with 19/26 socket tamping rod 40 cm. Long.

### Features:

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

Auto-zero on Reservoir Plain with Stop-cock