

Chromatography Columns With sintered Disc



Product Categories: [Chromatography Vials Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#)

Product Page:

<https://www.labappara.com/product/chromatography-columns-sintered-disc/>

Product Description

Chromatography Columns With sintered Disc - Laboratory Glassware Equipment

Chromatography Columns with integral sintered disc and PTFE key stopcock.

Features:

Long lasting performance

Durability

Lightweight

Accurate dimensions

Sleek design

Robust construction

Very high durability

Longer functionality

Auto-zero on Reservoir Plain with Stop-cock

Effective Length cm- 15 cm, 30 cm, 45 cm, 60 cm

Disc Dia mm - 10 mm, 20 mm, 30 mm